

description	RehaWash Systems GmbH		
Product:	Washing machine Wheelchair series 313	Washing machine Wheelchair series B 309	bed washing machine series B 904
			
Productie locatie:			
			
features			
suitable for	wheelchairs, commode chairs, rollators, AD matr.	hospital beds, AD matr., wheelchairs, rollators etc.	hospital beds, wheelchairs, rollators etc.
capacity	1 wheelchair	1 hospital bed in transport position, 1 wheelchair	1 hospital beds, 3 wheelchairs
cleaning programs	6: duration 8 /10/12 min. Adjustable	6: duration 8 /10/12 min. Adjustable	6: duration 8 /10/12 min. Adjustable
disinfection programs	2	2	2
custom wash programs possible	yes	yes	yes
document printing	yes	yes	yes
with drying	yes	yes	yes
available with Condense	yes	yes	yes
washing temperature	40-45 °C	40-45 °C	40-45 °C
mobile	yes	yes	yes
water consumption	35-40 liter	35-40 liter	35-40 liter
tape water	min. 2 bar	min. 2 bar	min. 2 bar
energy consumption	0,2 kW per wash cycle	0,2 kW per wash cycle	0,2 kW per wash cycle
energy supply	230 volt / 400 volt; 16A / 32A	230 volt / 400 volt; 16A / 32A	230 volt / 400 volt; 16A / 32A
operation			
Toutchtasten	yes	yes	yes
Display of information	yes	yes	yes
Technical ads selectable	yes	yes	yes
Counter for wash cycles	yes	yes	yes
Interface for information	yes	yes	yes
Roller door	by hand / automatic	by hand / automatic	by hand / automatic
Easy cleaning of the filter	yes	yes	yes
Filter material	stainless steal	stainless steal	stainless steal

description	RehaWash Systems GmbH		
	Product:	Washing machine Wheelchair series 313	Washing machine Wheelchair series B 309
security			
Emergency stop switch	yes	yes	yes
Automatischer stop when the door is open	yes	yes	yes
Display detergent empty	yes	yes	yes
Noise level in dB (A)	68	68	68
Disinfection test report	yes	yes	yes
dimensions			
dimensions l * w * h	1400 x 1640 x 2000 mm	1400 x 1640 x 2500 mm	2500 x 2050 x 2500 mm
dimensions l * w * h, inside	1350 x 860 x 1400 mm	1350 x 860 x 1900 mm	2450 x 1350 x 1900 mm
Dimensions selectable / adjust	yes	yes	yes
Delivery in modules / assembly at the customer	yes	yes	yes
lower it into the ground	yes	yes	yes
entrance height	min. 0 - max. 340 mm	min. 0 - max. 340 mm	min. 0 - max. 340 mm
maintenance			
life of at least	6-8 years	6-8 years	6-8 years
Number of wash cycles / year	3.000	3.000	3.000
Display of information	yes	yes	yes
Filter easy to clean - no tools	yes	yes	yes
Materials for filter	stainless steel	stainless steel	stainless steel
Software changes - at the customer's	yes	yes	yes
error message	yes, indicator light / horn	yes, indicator light / horn	Yes, indicator light / horn
Other auxiliary devices			
transport trolley	yes	yes	yes
A trolley according to customer requirements	yes	yes	yes
Normen			
CE, Medical Device class II a	93/42 EEC, ISO 15883	93/42 EEC, ISO 15883	93/42 EEC, ISO 15883
TÜV	yes	yes	yes
EMV	yes	yes	yes
Garantie			
year	1	1	1



RehaWash Systems GmbH
 Berliner Strasse 104
 02943 Weisswasser
 T +49 (0) 3576 223050
info@rehawash.com
www.rehawash.com