



Automatic washing line



AUTOMATIC WASHING LINE SPECIFICATIONS

Measurements

- Length: 5000
- Width: 2110
- Height: 2500

Capacity

- Bedstead: 1
- Bedstead foldable: 2
- Wheelchair: 3-4
- AD mattresses: 4
- Otherwise: 4

Operating parameters

- Operation voltage: 2 x 400 V
- Electrical power: 24 kWh
- Water pressure: 3-5 bar
- Water consumption: 35 - 40 l
- Sound level: 63 dB