



# Serie 313



## Additional info

A compact cleaning and disinfecting system for wheelchairs and rehabilitation accessories. The optional wheel washing device can be used to clean the large wheelchair rear wheels. The use of a special AD cleaning system enables the cleaning and disinfecting of anti-decubitus mattresses. Furthermore, optional cleaning baskets can be supplied for the treatment of smaller rehabilitation accessories.

# 313

## SPECIFICATIONS

### Measurements

- Length: 1400
- Width: 1640
- Height: 2000

### Capacity

- Bedstead: -
- Bedstead foldable: -
- Wheelchair: 1
- AD mattresses: 1
- Otherwise: 1
- Roller door: 1

### Operating parameters

- Operation voltage: 230 V / 400 V
- Electrical power: 2,2 kWh / 24 kWh
- Water pressure: 3-4 bar
- Water consumption: 35 - 40 l
- Sound level: 63 dB

## OPTIONS

### Accessories

- AD cleaning system: Yes
- Ramp: Yes
- Drip tray: Yes
- Pressure washer: Yes
- Washing device bike: Yes
- Cleaning basket: Yes
- Rail system: -
- Transport vehicle: -
- Drying system: -
- Water heating: Yes
- Additional heating (water): Yes

# 313 D

## SPECIFICATIONS

### Measurements

- Length: 1400
  - Width: 1640
  - Height: 2000
- 

### Capacity

- Bedstead: -
  - Bedstead foldable: -
  - Wheelchair: 1
  - AD mattresses: 1
  - Otherwise: 1
  - Roller door: 2
- 

### Operating parameters

- Operation voltage: 230 V / 400 V
  - Electrical power: 2,2 kWh / 24 kWh
  - Water pressure: 3-4 bar
  - Water consumption: 35 - 40 l
  - Sound level: 63 dB
- 

## OPTIONS

### Accessories

- AD cleaning system: Yes
  - Ramp: Yes
  - Drip tray: Yes
  - Pressure washer: No
  - Washing device bike: -
  - Cleaning basket: -
  - Rail system: -
  - Transport vehicle: -
  - Drying system: Yes
  - Water heating: -
  - Additional heating (water): -
-