



Serie 904 / B 904 D



Additional info

The RehaWash wash cell with integrated cleaning and disinfecting system is a validated and certified hygienic process for cleaning and disinfecting.

It clearly exceeds the requirements of the Medical Products Act for the safe treatment of reusable aids – demonstrated by hygienic analyzes based on a complex series of tests.

B 904

SPECIFICATIONS

Measurements

- Length: 2500
- Width: 2110
- Height: 2500

Capacity

- Bedstead: 1
- Bedstead foldable: 2
- Wheelchair: 3-4
- AD mattresses: 4
- Otherwise: 4

Operating parameters

- Operation voltage: 400 V
- Electrical power: 2,2 kWh / 24 kWh
- Water pressure: 3-5 bar
- Water consumption: 35 - 40 l
- Sound level: 63 dB

OPTIONS

Accessories

- AD cleaning system: Yes
- Ramp: Yes
- Drip tray: No
- Pressure washer: Yes
- Washing device bike: No
- Cleaning basket: Yes
- Rail system: Yes
- Transport vehicle: Yes
- Drying system: Yes
- Water heating: Yes
- Additional heating (water): Yes