



# Serie B 309 / B 309 D



## Additional info

A compact cleaning and disinfecting system for care and hospital beds in transport position (foldable beds). It is also possible to clean and disinfect wheelchairs, toilet chairs, anti-decubitus mattresses or rehab accessories in these wash cells.

# B 309

## SPECIFICATIONS

### Measurements

- Length: 1400
  - Width: 1640
  - Height: 2500
- 

### Capacity

- Bedstead: -
  - Bedstead foldable: 1
  - Wheelchair: 1
  - AD mattresses: 1
  - Otherwise: 1
  - Roller door: 1
- 

### Operating parameters

- Operation voltage: 230 V / 400 V
  - Electrical power: 2,2 kWh / 24 kWh
  - Water pressure: 2-3 bar
  - Water consumption: 35 - 40 l
  - Sound level: 63 dB
- 

## OPTIONS

### Accessories

- AD cleaning system: Yes
  - Ramp: Yes
  - Drip tray: No
  - Pressure washer: Yes
  - Washing device bike: No
  - Cleaning basket: Yes
  - Rail system: Yes
  - Transport vehicle: -
  - Drying system: Yes
  - Water heating: Yes
  - Additional heating (water): Yes
-

# B 309 D

## SPECIFICATIONS

### Measurements

- Length: 1400
- Width: 1640
- Height: 2500

### Capacity

- Bedstead: -
- Bedstead foldable: 1
- Wheelchair: 1
- AD mattresses: 1
- Otherwise: 1
- Roller door: 2

### Operating parameters

- Operation voltage: 230 V / 400 V
- Electrical power: 2,2 kWh / 24 kWh
- Water pressure: 3-4 bar
- Water consumption: 35 - 40 l
- Sound level: 63 dB

## OPTIONS

### Accessories

- AD cleaning system: Yes
- Ramp: Yes
- Drip tray: Yes
- Pressure washer: No
- Washing device bike: -
- Cleaning basket: -
- Rail system: -
- Transport vehicle: -
- Drying system: Yes
- Water heating: -
- Additional heating (water): -